

## 90146 Drop Down Rail with Adjustable Leg User Instructions



CE

The Deluxe Drop Down Grab Rail with adjustable leg is designed to provide support to those who require additional help when sitting or raising themselves from the toilet. The rail needs to be fitted to a wall adjacent to the toilet with fixings that are suitable to the substrate of the wall to which it is being attached.

### Features

- Stylish but strong wall fixing rail (with 5 fixing holes), designed to aid toileting
- Locks positively into the wall mounted cast aluminium alloy bracket, to provide compact unit
- Easy lift up/fold down mechanism
- Better grip comfort from oval profile hand rail
- Soft feel grey end cap
- Hard wearing nylon modified powder coated steel

## Specifications

- Length of Rail 70cm (27.5") (from wall to tip)
- Max Distance from Wall (when upright) 14.5cm (5.75")
- Mounting Bracket 26cm x 10cm (10.25" x 4")
- Adjustable Leg 69-78cm (27-30")
- Max user weight with leg 190Kg (30st)

## Cleaning Instructions

The rail should be cleaned with a anti-bacterial cloth or warm soapy water after every use & dried with a clean dry, soft cloth.