

Bed Grab Rail 90325



Description:

The 90325 Bed Grab Rail provides support when getting into and out of bed. It is designed for use on divan style beds and is suitable for use on single, double or kingsize beds. This all-metal grab rail is simple to assemble and dismantle without the need for any tools. It can be flattened for ease of storage and is a handy device to take with you on your holidays.

Specification:

Width	365 mm (14.5")
Length	620 mm (24.5")
Handle Height	470 mm (18.5")
Product Weight	3.1 kg
Max User Weight	160 kg (25 stone)

Cleaning:

The Integrity 90325 Bed Grab rail can be cleaned using a damp cloth with any mild detergent. **Do Not** use any type of abrasive on this product. After cleaning, dry with a clean dry cloth.