


90397

***Postural Wedge -
Short Instruction***



Illustrations

Illustration 1



Materials

Outer material: Polyurethane-coated nylon, water-resistant.
Postural Wedge: Soft & Stable Cold Foam



Cleaning

The product should be checked regularly, preferably each time it is used, and especially after washing.

Check that there is no damage to seams or fabric.



Never try to repair a product yourself. Never use a defective product.

Washing instructions:

Remove foam from cover before washing.

Do not add fabric conditioner.

Can also be disinfected using LT-CIN (Isopropyl Alcohol 72%).

This product can be cleaned with a solvent-free cleaning agent with a pH level between 5 - 9 or a 70% disinfection solution.

Key to Symbols:



Washing Temp.



DO NOT Iron



Never leave on the floor



Tumbledry



Handwash



Do not dryclean



Dispose of Responsibly



Wipe off



Do not bleach



Read the Instructions



Risk of sliding down