

Telescopic Suction Grab Bar Part Numbers: 90631, 90632 & 90633



Description

These bars can be installed anywhere for your safety without having to drill unsightly holes in your wall tiles or walls. It installs in seconds without tools. Movable at leisure and ideal for your travels and trips. Press the two side levers and it holds very firmly on all smooth surfaces (tiles, earthenware, glass surface, lacquered paint, PVC surface of shower cubicles...)

- Bar adjustable over 12 cm
- Very strong rubber suction cups
- Secure suction cups by lowering the lever
- Maximum lifting weight: 40 kg

3 different lengths :

Length	Part Number
44-56 cm	90631
55-68 cm	90632
65-78 cm	90633

- Suction cup diameter: 12 cm